

Polybutadiene (or rubber) is added during polymerisation of styrene to increase the flexibility of polystyrene polymers. It is important to closely monitor the amount of polybutadiene added, so that the flexibility of the polystyrene polymer may be regulated. NMR can provide an effective and rapid method for the regulation of polybutadiene content in Impact Polystyrene.

### Method

NMR has a number of advantages over other techniques:

- NMR is very stable over the long term and rarely needs calibration adjustment.
- The measurement time is short (typically 32 seconds), allowing rapid sample throughput.
- The NMR technique is non-destructive, so the same sample may be measured several times before being analysed by other techniques.
- NMR is insensitive to air voids between polystyrene granules.
- Measurement precision is good compared to wet chemical methods.
- NMR does not require the use of hazardous solvents.
- Both weighing and non-weighing methods are available for this application

### Calibration and Results

12 samples of impact polystyrene were used for the calibration of polybutadiene content. Each sample was weighed into a tared 26mm glass tube before measurement at room temperature in the **MQC+**. The resulting calibration graph is shown in Figure 1.




Table 1 shows that the repeatability (or precision) of the NMR data is excellent.

Repeat number	Polybutadiene Content (%)
1	6.79
2	6.81
3	6.80
4	6.81
5	6.79
6	6.78
7	6.81
8	6.79
9	6.79
10	6.80
Average	6.80
Standard Deviation	0.01

**Table 1:** NMR measurement repeatability.

**Figure 1:** Calibration graph of NMR signal per unit mass against polybutadiene content obtained by a reference method. Correlation Coefficient  $r = 1$ , Standard Deviation = 0.15.



The Business of Science®

## Recommended Instrument

The **MQC+** with a 0.55 Tesla (23 MHz hydrogen resonance) magnet, fitted with a 26mm diameter (14ml sample) probe is a suitable instrument for this application. The Rubber in Impact Polystyrene package consists of:

- **MQC+** with a built-in computer operating Microsoft® Windows®10 (no separate PC is required).
- **MultiQuant** software including **RI Calibration**, **RI Analysis**, and the **EasyCal** 'Rubber in Polystyrene' application.
- 26mm diameter glass tubes.
- User Manuals.

**In addition to this package you may also require:**

- A precision balance.

**The instrument offers advantages over others on the market:**

- High signal sensitivity.
- Small benchtop footprint.
- Low maintenance.
- Recyclable sample tubes, lowering consumable costs.
- Minimal sample preparation.

**Note:** Other instruments/packages are available for the analysis of larger or smaller quantities of sample. Please contact Oxford Instruments for details.



visit [www.oxford-instruments.com/mqc](http://www.oxford-instruments.com/mqc) for more information or email: [magres@oxinst.com](mailto:magres@oxinst.com)

This publication is the copyright of Oxford Instruments and provides outline information only which (unless agreed by the company in writing) may not be used, applied or reproduced for any purpose or form part of any order or contract or be regarded as a representation relating to the products or services concerned. Oxford Instruments' policy is one of continued improvement. The company reserves the right to alter, without notice, the specification, design or conditions of supply of any product or service. Oxford Instruments acknowledges all trademarks and registrations. © Oxford Instruments plc, 2018. All rights reserved. Part no. PIP-03-18



*The Business of Science®*

